

Assay Protocol for Endogenous GPCR Activation in Human Pulmonary Artery Smooth Muscle Cells (PASMCs) Expressing AEQ-GFP in Mitochondria (cAP-0101AEQ-GFP-Mito)

Purpose

This protocol describes a **functional bioluminescence assay** to measure endogenous GPCR activation in **human pulmonary artery smooth muscle cells (PASMCs) stably expressing mitochondrial-targeted AEQ-GFP (cAP-0101AEQ-GFP-Mito)** using **aequorin-based calcium mobilization**.

Materials Required

1. **Cells:** Human PASMCs expressing AEQ-GFP in mitochondria (cAP-0101AEQ-GFP-Mito)
 2. **Culture Medium:** SMCs Growth Medium (cAP-24)
 3. **GPCR Ligands:** Based on commonly expressed GPCRs in PASMCs (see list below)
 4. **Coelenterazine h (or native coelenterazine):** For reconstitution of aequorin (5 μ M final concentration)
 5. **Assay Buffer:** HBSS + 0.1% BSA (pH 7.4)
 6. **CaCl₂ (for calcium reintroduction if needed)**
 7. **Luminometer (e.g., Berthold MicroLumat LB 96P or equivalent)**
 8. **White 96-well or 384-well microplates (luminescence-compatible)**
 9. **Trypsin/EDTA (for cell detachment)**
 10. **Probenecid (optional, 5 mM, to prevent dye leakage)**
-

Procedure

Day 1: Cell Seeding

1. **Thaw and Culture Cells:**
 - Recover cells from liquid nitrogen and culture in **SMCs Growth Medium (cAP-24)** at **37°C, 5% CO₂** until ~80% confluent.
 - Passage cells using **trypsin/EDTA** as needed.
2. **Seed Cells for Assay:**

- Detach cells and seed in a **white 96-well plate** at **50,000–80,000 cells/well** in **100 μ L growth medium**.
 - Allow cells to adhere overnight (~16–24 hr).
-

Day 2: Aequorin Reconstitution & GPCR Stimulation

- 1. Prepare Assay Buffer:**
 - **HBSS + 0.1% BSA** (pH 7.4).
 - (Optional) Add **5 mM probenecid** to prevent dye leakage.
 - 2. Reconstitute Aequorin with Coelenterazine h:**
 - Dilute **coelenterazine h** in assay buffer to **5 μ M** (final concentration).
 - Remove growth medium and add **100 μ L/well of coelenterazine h solution**.
 - Incubate **2–4 hr at 37°C** (or **overnight at RT** for better signal stability).
 - 3. Prepare GPCR Ligands:**
 - Dilute ligands in assay buffer at **10X final desired concentration** (accounting for 1:10 dilution upon addition to cells).
 - **Example concentrations for PASMIC-expressed GPCRs:**
 - **Endothelin-1 (ETAR/ETBR):** 10–100 nM
 - **Angiotensin II (AT1R):** 100 nM–1 μ M
 - **Thrombin (PAR1/PAR4):** 0.1–1 U/mL
 - **Serotonin (5-HT_{2B} receptor):** 1–10 μ M
 - **ATP (P2Y receptors):** 1–10 μ M
 - **Prostaglandin F₂ α (FP receptor):** 100 nM–1 μ M
 - **Sphingosine-1-phosphate (S1PR1–3):** 100 nM–1 μ M
 - **Urotensin-II (UT receptor):** 10–100 nM
 - 4. Run the Assay (Luminometer Setup):**
 - Pre-equilibrate plate at **37°C** for 10 min.
 - Set up luminometer to inject **10 μ L of 10X ligand** per well (total volume = 110 μ L).
 - Measure luminescence **immediately after ligand addition** (1–2 sec intervals for 30–60 sec).
-

Data Analysis

- **Peak Luminescence:** Measure maximum signal after ligand addition.
- **Area Under the Curve (AUC):** Integrate signal over time for total response.

- **Normalization:** If needed, normalize to baseline or positive control (e.g., ionomycin/Ca²⁺ ionophore for maximum response).
-

Commonly Expressed GPCRs in Human PSMCs

PASMCs express a variety of endogenous GPCRs, including:

- **Endothelin receptors (ETAR, ETBR)** – Vasoconstriction (ET-1)
 - **Angiotensin II receptor (AT1R)** – Vasoconstriction (Ang II)
 - **Protease-activated receptors (PAR1, PAR4)** – Proliferation (thrombin)
 - **Serotonin receptors (5-HT2B, 5-HT1B)** – Vasoconstriction (5-HT)
 - **Purinergic receptors (P2Y1, P2Y2, P2Y6, P2Y12)** – ATP/ADP-mediated signaling
 - **Prostanoid receptors (FP, TP, EP1–4)** – PGF2 α /TXA2/PGE2 signaling
 - **Sphingosine-1-phosphate receptors (S1PR1–3)** – S1P-mediated migration/proliferation
 - **Urotensin-II receptor (UT)** – Potent vasoconstrictor (U-II)
 - **Muscarinic receptors (M3, M5)** – Acetylcholine-mediated contraction
-

Troubleshooting

- **Low Signal:** Increase coelenterazine incubation time or concentration.
 - **High Baseline Noise:** Reduce coelenterazine exposure to light (light-sensitive).
 - **No Response:** Verify GPCR expression in PSMCs; test positive control (e.g., ATP or ionomycin).
-

Conclusion

This protocol allows functional assessment of **endogenous GPCR activation** in **PASMCs (cAP-0101AEQ-GFP-Mito)** via **mitochondrial Ca²⁺-dependent aequorin luminescence**.