

Penton-Buford's Mouse Endothelial Growth Medium
ORDERING INFORMATION

Product Name	Penton-Buford's Mouse Endothelial Growth Medium (PBMEGM)
Catalog No:	cAP-64
Size:	500ml
Storage:	4°C for 2 months

Product Description

Penton-Buford's Mouse Endothelial Growth Medium (PBMEGM) is specifically formulated for long term culture of mouse primary endothelial cells.

Angio-Proteomie has exclusive license right for the sales and distribution of the OBMEGM worldwide.

Mouse primary ECs can be passaged and cultured for up to 20 passages when cultured using PBMEGM. It has been tested also works for certain Rat and human ECs.


Shipping Condition: Ambient temperature or blue ice (seasonally)

Storage Condition: PBMEGM should be stored at 4°C. A change in color or appearance of precipitate may indicate deterioration.


Shelf Life: 2 months from the date of receipt under proper storage condition.

1. Ready to use.
2. Avoid frequent temperature changes to the entire bottle of medium.
3. We recommend medium changing every 2 days for regular cell expansion culture.

PMVEC -MOUSE



PMVEC - RAT


Related Products:

Quick Coating Solution	cAP-01	240ml	Angio-Proteomie
Cell Freezing Solution (FBS)	cAP-22	50ml	Angio-Proteomie
Cell Freezing Solution (Non-FBS)	cAP-22B	50ml	Angio-Proteomie
HBSS w/o Ca ²⁺ , Mg ²⁺	cAP-11	100ml	Angio-Proteomie
Trypsin/EDTA Solution	cAP-23	100ml	Angio-Proteomie
Trypsin Neutralization Solution	cAP-28	100ml	Angio-Proteomie
ITS (100x)	cAP-26	10ml	Angio-Proteomie
L-Glutamine-MAXIMUM (100x)	cAP-27	100ml	Angio-Proteomie
Human Plasma Fibronectin Solution	cAP-42	1mg/ml	Angio-Proteomie
Bovine Type I Collagen Solution	cAP-17	100mg	Angio-Proteomie

THESE PRODUCTS ARE FOR RESEARCH USE ONLY

Caution: Handling human and animal tissue derived products is potentially bio-hazardous. Although each cell strain is tested negative for HIV, HBV and HCV DNA, or pathogens, diagnostic tests are not necessarily 100% accurate; therefore proper precautions must be taken to avoid inadvertent exposure. Always wear gloves and safety glasses when working with these materials. Never mouth pipette. We recommend following the universal procedures for handling products of human origin as the minimum precaution against contamination.