

Keratinocyte Growth Medium
ORDERING INFORMATION

Product Name	Keratinocyte Growth Medium
Catalog No:	cAP-67
Size:	500ml
Storage:	4°C

Product Description

Keratinocyte Growth Medium (KeraGM) is freshly prepared by mixing the following two components: (1) Keratinocyte Basal Medium and Epithelial Growth Supplement (Not Available Separately), which contains the essential nutrients for Keratinocyte growth (containing recombinant growth factors, and 1x penicillin and streptomycin).

Application: This medium is special designed for Keratinocyte cell expansion.

Shipping Condition: Ambient temperature or blue ice (seasonally)

Storage Condition: The KeraGM should be stored at 4°C. A change in color or appearance of precipitate may indicate deterioration.

Shelf Life: at least 1 months from the date of receipt under proper storage condition.

Related Products:

Quick Coating Solution	cAP-01	240ml	Angio-Proteomie
Cell Freezing Solution (FBS)	cAP-22	50ml	Angio-Proteomie
Cell Freezing Solution (Non-FBS)	cAP-22B	50ml	Angio-Proteomie
HBSS w/o Ca ²⁺ , Mg ²⁺	cAP-11	100ml	Angio-Proteomie
Trypsin/EDTA Solution	cAP-23	100ml	Angio-Proteomie
Trypsin Neutralization Solution	cAP-28	100ml	Angio-Proteomie
ITS (100x)	cAP-26	10ml	Angio-Proteomie
L-Glutamine-MAXIMUM (100x)	cAP-27	100ml	Angio-Proteomie
Human Plasma Fibronectin Solution	cAP-42	1mg/ml	Angio-Proteomie
Bovine Type I Collagen Solution	cAP-17	100mg	Angio-Proteomie
Doxycycline Solution (1000 x)	cAP-45	10ml	Angio-Proteomie

THESE PRODUCTS ARE FOR RESEARCH USE ONLY

Caution: Handling human and animal tissue derived products is potentially bio-hazardous. Although each cell strain is tested negative for HIV, HBV and HCV DNA, or pathogens, diagnostic tests are not necessarily 100% accurate; therefore proper precautions must be taken to avoid inadvertent exposure. Always wear gloves and safety glasses when working with these materials. Never mouth pipette. We recommend following the universal procedures for handling products of human origin as the minimum precaution against contamination.